

# SL19A    Advanced Slit Lamp Microscope



## Features

- New developed parallel optical path  
Microscpe provides higher depth of view  
and larger-angle observation
- All the lens are multi-coated, mildewproof  
and anti-reflection
- LED and halogen light soures are optional

## Specifications

Model	SL19A
MICROSCOPE	
Microscope type	Parallel stereomicroscope
Eyepiece	12.5X
Magnification	6X, 10X, 16X, 25X, 40X
Field of View	φ5.5, φ9, φ14 ,φ23, φ35(mm)
Adjustable Pupil Distance	55-72mm
Adjustable Diopter	> ±7D
SLIT	
Slit Width	0-14mm adjustable
Slit Length	1-14mm adjustable
Diameter of Aperture	φ0.2, φ2, φ3, φ5, φ8, φ10, φ14
Slit Angle	0-180°
Slit Inclination	5°, 10°, 15°, 20°
Filter	Heat absorption, Grey, Redfree, Cobalt blue
Lighting bulb	3.3V 10W LED Lamp or halogen lamp
Fixed view lamp	Red LED
POWER	
Input voltage	AC110V/60Hz or AC220V/50Hz
Input power	80VA
Chin-rest Vertical Movement	≥80mm
Base Longitudinal Movement	≥90mm
Base Lateral Movement	≥100mm
Base Vertical Movement	≥30mm
Joystick	with shutter button