

SL1X2 (2steps) Slit lamp microscope



Features

- Simple slit lamp microscope with 2 steps magnifications.
- Excellent optical materials gives stable performance.
- LED Illumination.

Specifications

Model	SL1X2
MICROSCOPE	
Microscope type	Binocular converging stereomicroscope
Eyepiece	10X
Magnification	10X, 16X
Field of View	φ14, φ17.5 (mm)
Adjustable Pupil Distance	55-72mm
Adjustable Diopter	> ±5D
SLIT	
Slit Width	0~14mm adjustable
Slit Length	1-14mm adjustable
Spot diameter(mm)	φ0.3, φ2, φ5.5, φ9, φ14
Slit Angle	0-180°
Filter	Heat absorption, Grey, Redfree, Cobalt blue
Lighting bulb	3.3V,10W LED lamp
Fixed view lamp	Red LED
POWER	
Input voltage	AC110V/60Hz or AC220V/50Hz
Input power	80VA
Optional Accessories	Instrument table, Applanation tonometer, Smartphone adaptor