

PROBABLY THE FASTEST OBJECTIVE VISION SCREENER


- ⦿ Objective and accurate refractive assessment
- ⦿ Ideal for challenging patients such as infants, disabled patient
- ⦿ Revolutionary design brings superior user experience

EASY OPERATION THAT YOU CAN MASTER IN SECONDS



Step 1

Simply put the screener against the forehead



Step 2

Centre the cross on the pupil



Step 3

Save, print or share the report

Sphere range	-10.0D ~ +10.0D
Sphere tolerance	±0.5D
Cylinder range	-5.0D ~ +5.0D
Cylinder tolerance	±0.5D
Axis of astigmatism range	1° ~ 180°
Data acquisition time	1s/eye
External	External must connect to PC and IPAD
Net weight	360g

Specification subject to change without notice